

Change Log
– ezColony®

Version 5.0.1		<i>Release Date:</i>
Fix	Cage view fails to load if preferred Animal Line deleted	

Version 5.0.0		<i>Release Date:</i>
-	New release	
-	New streamlined UI	
-	Expanded Advanced Search	

Version 4.1.17		<i>Release Date:</i>
#2735	Feature: Ability to see lost animal location in archive record (Novonordisk request)	

Version 4.1.16		<i>Release Date:</i>
Add	Login Message	
Add	Single – so option	

Version 4.1.15		<i>Release Date:</i>
Fix	Improved LDAP integration	
Add	Changed inactivity period	

Version 4.1.14		<i>Release Date:</i>
Add	New custom icons	

Version 4.1.12		<i>Release Date:</i>
Fix	Quick search. Sometimes it may return extra records	
Fix	API search for archived animals	

Version 4.1.11		<i>Release Date:</i>
Add	Expanded 'animals' API call	

Version 4.1.10		<i>Release Date:</i>
Fix	RadioCheck box control	

Version 4.1.9		<i>Release Date:</i>
Add	Internal ID to the animal report	

Version 4.1.8		Release Date:
Fix	Save object dialogs that in some cases show incorrect error message (Chrome 31)	

Version 4.1.7		Release Date:
Bug		
#1258	ezColony®, 4.1.4RC - VM Console - Sync with LDAP - post execution, users created with Internal Authentication are not disabled	
#1355	ezColony® 4.1.6 Checkbox User-Defined field doesn't allow unchecking	

Version 4.1.5		Release Date:
Fix	Fixed mating procedure for MS SQL DB	

Version 4.1.4		Release Date:
Fix		
#1229	ezColony® - Selecting Animals - Individual Checkboxes Do Not Work	
#1242	ezColony®, 4.1.4RC - Notch import does not work when Inverted Orientation is selected	
#1248	ezColony® 4.1.4 - User Preferences - Changing The Color Theme	
#1249	ezColony® 4.1.4. - Editing User Groups - Users List Is Emptied When Editing	
#1251	ezColony® 4.1.4 - Cage Types - Editing Cage Changes Created Date	
#1252	ezColony® 4.1.4. - IE8 UDFs - Windows With Many UDFs Don't Scroll	
#1255	ezColony®, 4.1.4RC - VM Console - attempts to change Display Name Setting did not succeed	
Add		
#1137	ezColony® 4 - Hiding Specific Fields - Global Effect	
#1223	To add ability to delete a virtual animal(s) from a Pick List	

Version 4.1.3		<i>Release Date:</i>
Fix		
#783	ezColony® 3.5.10 - Chinese - Add New User - No Role Dropdown	
#1039	ezColony® 4.1.0 - Advanced Search – Sorting	
#1040	ezColony® 4.1.0 - Advanced Search Archive – Sorting	
#1201	Wrong animals count for users who don't have 'System Administration' right	
#1206	ezColony®, 4.1.3 - import animals - Date In (Date of Birth) does not import correctly	
#1210	ezColony®, 4.1.3 - Quick Search tool - picks up records without the search term	
-	Quick search. Sometimes it may return extra records	
-	API search for archived animals	
Add		
#1196	Implement Advanced Alerts similar to FreezerPro®	
-	Add Custom Rack coordinates option	

Version 4.1.2		<i>Release Date:</i>
Fix	Animal Line Creation	

Version 4.1.1		<i>Release Date:</i>
Add	Custom (NIH) UL Number report	

Version 4.1.0		<i>Release Date:</i>
Fix		
	Group edit for date and list UDFs	
	Auto-increment UDF	
	Animals, Cages batch search on UDFs	
Add		
	Litters UDFs	
	Support for IE10	
	Advance Search on Archived Animals	
	New custom cage printing profile	
	'None' ear notch schema	

Version 4.0.1		<i>Release Date:</i>
Fix	Alerts on MS SQL	
Fix	Animal view in case of a cage having an ' symbol in the name	

Version 4.0.0		<i>Release Date:</i>
Fix	CSV export opens in new browser tab	
Fix	Creating empty user groups	
Add	New alert for mating cages with animals older than 9 months	

Version 3.5.7		<i>Release Date:</i>
Fix	Backup procedure	

Version 3.5.6		<i>Release Date:</i>
Fix	Online Activation	

Version 3.5.5		<i>Release Date:</i>
Fix	Import files with bad dates	

Version 3.5.4		<i>Release Date:</i>
Fix	File UDF	
Fix	CSV reports	
-	UI improvements	

Version 3.5.2		<i>Release Date:</i>
Fix		
-	Advanced Search – in rare cases, choice type UDF values might not refresh in Query Builder	
-	Actions menu on iPad	
-	Editing when User type UDF is removed	
-	Admins can Edit Animals in Archive	
-	Admins can move Animals from Archive	
-	Admins can change disposition date in Archive	
New	Option to move animals by scanning cage barcodes	

Version 3.5.1		<i>Release Date:</i>
-	Improved LDAP authentication code	
	Removed unused Cell Line object field	

Version 3.5.0		<i>Release Date:</i>
Fix	Full IE 9.0 support	
New	Ear Notch Schema	

Version 3.2.5		<i>Release Date:</i>
Fix	IE 9.0	

Version 3.2.3		<i>Release Date:</i>
-	Option to Import Genotypes and Ear Notches added	

Version 3.2.2		<i>Release Date:</i>
-	New update procedure	

Version 3.2.1		<i>Release Date:</i>
Fix	Bug with cell line list not showing	
Fix	Bug with deleting UDFs	
Add	Type-a-head for cell lines and users	

Version 3.2.0		<i>Release Date:</i>
New	New UDF: ezColony® User	
New	New UDF: ezColony® Animal Line	
-	Support for harem breeding	
-	Remember last mating cage for female animals	
-	Option to rearrange UDFs in entry forms	
-	Option to rotate animal ear-notch images	

Version 3.1.4		<i>Release Date:</i>
Fix	Fix for mice editing for Oracle DB	
Add	Scrolling for all dialog with UDF	
-	Inactivity dialog fixed	

Version 3.1.3		<i>Release Date:</i>
-	Renamed columns that had Oracle keyword names	

Version 3.1.2		<i>Release Date:</i>
Fix	Removing animals from the cage	
Fix	Time display in grids	

Version 3.1.0		<i>Release Date:</i>
Add	Mouse operation icons into the cage list view	
Add	Batch Search of animals and cages	
-	Animal Pick Lists (like in FreezerPro®)	
-	Show cage barcode in animals lists	
-	Show cage barcode in the rack view	
-	Option to add User Defined Fields to Requests (Enterprise)	
-	Send requests to User Groups (not individual users)	
-	Option to assign access permissions to UDFs (Enterprise)	
-	3 new cage type icons	
-	Bug fixes: Proper message when restriction is in effect (moving/adding animals to cages)	

Version 3.0.12		<i>Release Date:</i>
-	Rack height extension to 12 (from 10)	

Version 3.0.11		<i>Release Date:</i>
Fix	nameservers setting in Network configuration dialog	
Fix	Default role for users created from LDAP	

Version 3.0.10		<i>Release Date:</i>
-	Print Cage cards from rack right-click menu	
-	Fix for reports in archive (when run from animal line)	

Version 3.0.9		<i>Release Date:</i>
Fix	Delete an ACL record if an object it belongs to is destroyed	
Fix	Respect ACLs in cage, litter lists	

Version 3.0.7		<i>Release Date:</i>
Fix		
-	For saving Animal Lines (duplicated genotypes)	
-	For All Animals - show all live animals	
-	When clicking the sex icon the menu appears in the top left corner	
-	Sorting by animal line not always works	

Version 3.0.6		<i>Release Date:</i>
-	Cage number in search fixed	
-	Cage printed fixed	

Version 3.0.5		<i>Release Date:</i>
-	NEW ezColony® 2011 Release	